
Novuqare[®] SealMateP40/P20

High-quality validatable rotary sealers. With touchscreen and optional printer.



Novuqare is making waves in the sealer market with the new SealMate range of products which are:

- High quality
- Extremely reliable
- User friendly
- Easy to service

The SealMate is the result of Novuqare's long experience in designing and building premier rotary sealers. Sealers from Novuqare have been engineered and tested to meet the demanding conditions of hospital sterilisation departments. Ten years of experience and continuous technology development guarantee very reliable and easy-to-service products.

With the arrival of Novuqare's SealMate heat sealers, a new benchmark in quality sealers has been set.

	Seal Mate P40	Seal Mate P20
Validatable (ISO 11607-2)	\checkmark	1
DIN EN 868-5, DIN 58953-7	\checkmark	1
9 dot matrix printer	\checkmark	×
Touchscreen	\checkmark	1
Temperature adjustable	\checkmark	1
Speed adjustable	\checkmark	1
Pressure adjustable	\checkmark	1
Auto stand-by	\checkmark	1
Auto start	\checkmark	1
Slow start	\checkmark	1
Adjustable pouch insert depth	 Image: A set of the set of the	1
Electrical power	110-230V / 50-60 Hz / 550W	
Weight	10kg / 22lbs	
Dimensions (WxDxH)	580x255x133mm	

Novuqare

For more information please contact us: www.novuqare.com +31 (0) 73 782 01 10 sealers@novuqare.com



Innovations for healthier lives.